



## SAFETY DATA SHEET

### **SECTION 1 – PRODUCT AND COMPANY IDENTIFICATION**

- |            |   |   |
|------------|---|---|
| <b>1.1</b> | PRODUCT NAME:<br>CATALOG #:             | RIPA Buffer<br>BP-115   |
| <b>1.2</b> | RECOMMENDED USE:                        | Research Use or Further Manufacturing                           |
| <b>1.3</b> | COMPANY IDENTIFICATION:                 | Boston BioProducts, Inc.<br>425 Fortune Blvd. Milford, MA 01757 |
| <b>1.4</b> | EMERGENCY CONTACT:<br>TOLL FREE NUMBER: | 1-800-535-5053<br>1-888-881-8691                                |

### **SECTION 2 - HAZARDS IDENTIFICATION**

**2.1 Classification of the substance or mixture**  
**GHS Classification in accordance with 29 CFR 1910 (OSHA HCS)**

|      |  |            |
|------|--|------------|
| H228 | Flammable solids                                 | Category 2 |
| H302 | Acute toxicity, Oral                             | Category 4 |
| H332 | Acute toxicity, Inhalation                       | Category 4 |
| H315 | Skin irritation                                  | Category 2 |
| H318 | Serious eye damage                               | Category 1 |
| H335 | Specific target organ toxicity – single exposure | Category 3 |
| H401 | Short-term (acute) aquatic hazard                | Category 2 |
| H412 | Long-term (chronic) aquatic hazard               | Category 3 |

**Signal word:** Danger

**Hazard statement(s):**

|             |  |
|-------------|--|
| H228        | Flammable solid.                                   |
| H302 + H332 | Harmful if swallowed or inhaled.                   |
| H315        | Causes skin irritation.                            |
| H318        | Causes serious eye damage.                         |
| H335        | May cause respiratory irritation.                  |
| H401        | Toxic to aquatic life.                             |
| H412        | Harmful to aquatic life with long lasting effects. |

**Precautionary statements:**

|      |   |
|------|---|
| P210 | Keep away from heat/ sparks/ open flames/ hot surfaces. No smoking. |
| P240 | Ground/ bond container and receiving equipment.                     |
| P241 | Use explosion-proof electrical/ ventilating/ lighting/ equipment.   |
| P261 | Avoid breathing dust.   |



|                           |   |
|---------------------------|---|
| P264                      | Wash skin thoroughly after handling.  |
| P270                      | Do not eat, drink or smoke when using this product.   |
| P271                      | Use only outdoors or in a well-ventilated area.   |
| P273                      | Avoid release to the environment.   |
| P280                      | Wear protective gloves/ eye protection/ face protection.  |
| P301 + P312 + P330        | IF SWALLOWED: Call a POISON CENTER/ doctor if you feel unwell. Rinse mouth.   |
| P302 + P352               | IF ON SKIN: Wash with plenty of soap and water.   |
| P304 + P340 + P312        | IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER/ doctor if you feel unwell.   |
| P305 + P351 + P338 + P310 | IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POSISON CENTER/ doctor. |
| P332 + P313               | If skin irritation occurs: Get medical advice/ attention.   |
| P362                      | Take off contaminated clothing and wash before reuse.   |
| P370 + P378               | In case of fire: Use dry sand, dry chemical or alcohol-resistant foam to extinguish.  |
| P403 + P223               | Store in a well-ventilated place. Keep container tightly closed.  |
| P405                      | Store locked up.  |
| P501                      | Dispose of contents/ container to an approved waste disposal plant.   |

**SECTION 3 - COMPOSITION/INFORMATION ON INGREDIENTS**

**Composition:** Mixture of 50 mM Tris Base and Tris-HCl, 150 mM Sodium Chloride, 0.5% Sodium Deoxycholate, 0.1% SDS, 1% Nonidet P-40 Substitute, pH 7.40 ± 0.15.

| <b>Component</b>        | <b>CAS Number</b> | <b>Concentration (%)</b> |
|-------------------------|-------------------|--------------------------|
| Tris Base               | 77-86-1           | 0.09%                    |
| Tris-HCl                | 1185-53-1         | 0.66%                    |
| Sodium Chloride         | 7647-14-5         | 0.88%                    |
| Sodium Deoxycholate     | 302-95-4          | 0.50%                    |
| Sodium Dodecyl Sulfate  | 151-21-3          | 0.10%                    |
| Nonidet P-40 Substitute | 9016-45-9         | 1.00%                    |
| Water                   | 7732-18-5         | rest                     |

**SECTION 4 - FIRST-AID MEASURES**

**4.1 Description of first aid measures**

**General advice:** Consult a doctor and show this safety data sheet.

- i. **If inhaled:** Remove to fresh air and monitor breathing. If breathing becomes difficult, give oxygen. If breathing stops, give artificial respiration. Consult a doctor.
- ii. **In case of skin contact:** Immediately wash skin with copious amounts of soap and water for at least 15 minutes. Remove contaminated clothing and shoes and wash before reuse. Consult a doctor.



- iii. **In case of eye contact:** Flush with copious amounts of water for at least 15 minutes. Consult a doctor.
- iv. **If swallowed:** Rinse mouth with water. Do not induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Consult a doctor.

**4.2 Most important symptoms and effects, both acute and delayed:** To the best of our knowledge, the chemical, physical and toxicological properties have not been thoroughly investigated.

**4.3 Indication of immediate medical attention and special treatment needed:** Show this safety data sheet to the doctor in attendance. Immediate medical attention is required.

## **SECTION 5 - FIRE FIGHTING MEASURES**

### **5.1 Extinguishing media**

**Suitable extinguishing media:** Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide

**5.2 Special hazards arising from the substance or mixture:** No data available.

**5.3 Precautions for fire-fighters:** Wear suitable protective clothing to prevent contact with skin and eyes and self-contained breathing apparatus.

**5.4 Further information:** No data available

## **SECTION 6 - ACCIDENTAL RELEASE MEASURES**

**6.1 Personal precautions, protective equipment and emergency procedures:** Do not take action without suitable protective clothing - see section 8 of SDS. Evacuate personnel to safe areas. Ensure adequate ventilation. Avoid breathing vapors, mist, dust or gas.

**6.2 Environmental precautions:** Do not let product enter drains.

**6.3 Methods and materials for containment and cleaning up:** Cover spillage with suitable absorbent material. Using non-spark tools, sweep up material and place in an appropriate container. Decontaminate spill site with 10% caustic solution and ventilate area until after disposal is complete. Hold all material for appropriate disposal as described under section 13 of SDS.

**6.4 Reference to other sections:** For required PPE see section 8. For disposal see section 13.

## **SECTION 7 - HANDLING AND STORAGE**

**7.1 Precautions for safe handling:** Use in a chemical fume hood, with air supplied by an independent system. Avoid inhalation, contact with eyes, skin and clothing. Avoid the formation of dust and aerosols. Use in a well-ventilated area. Keep away from sources of ignition. Avoid prolonged or repeated exposure.

**7.2 Conditions for safe storage, including any incompatibilities:** Store in cool, well-ventilated area. Keep away from direct sunlight. Keep container tightly sealed until ready for use. Keep in a dry place. Store at 2 - 8 °C.

**7.3 Specific end use(s):** Use in a laboratory fume hood where possible.

## **SECTION 8 - EXPOSURE CONTROLS/PERSONAL PROTECTION**

### **8.1 Control parameters**

**Components with workplace control parameters:** Contains no substances with occupational exposure limit values.

### **8.2 Exposure Controls**

- i. **Appropriate engineering controls:** Use in a fume hood where applicable. Ensure all engineering measures described under section 7 of SDS are in place. Ensure laboratory is equipped with a safety shower and eye wash station.



**8.3 Personal protective equipment**

- i. **Eye/face protection:** Face shield and safety glasses Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).
- ii. **Skin protection:** Wear appropriate protective eyeglasses or chemical safety goggles as described by OSHA's eye and face protection regulations in 29 CFR 1910.133 or European Standard EN166
- iii. **Body Protection:** Wear appropriate protective clothing. Complete suit protecting against chemicals, The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.
- iv. **Respiratory Protection:** Follow the OSHA respirator regulations found in 29 CFR 1910.134 or European Standard EN 149. Use a NIOSH/MSHA or European Standard EN149 approved respirator if exposure limits are exceeded or if irritation or other symptoms are experienced.
- v. **Control of environmental exposure:** Do not let product enter drains

**SECTION 9 - PHYSICAL AND CHEMICAL PROPERTIES**

**Information on basic physical and chemical properties**

|   |                         |                                  |                   |
|---|-------------------------|----------------------------------|-------------------|
| <b>Appearance</b>                                     | Clear, Colorless Liquid | <b>Vapor Pressure</b>            | No Data Available |
| <b>Odor</b>   | No Data Available       | <b>Vapor Density</b>             | No Data Available |
| <b>Odor Threshold</b>                                 | No Data Available       | <b>Relative Density</b>          | No Data Available |
| <b>pH</b>   | 7.40                    | <b>Water Solubility</b>          | No Data Available |
| <b>Melting / Freezing Point</b>                       | No Data Available       | <b>Partition Coefficient</b>     | No Data Available |
| <b>Initial Boiling Point Range</b>                    | No Data Available       | <b>Auto-Ignition Temperature</b> | No Data Available |
| <b>Flash Point</b>                                    | No Data Available       | <b>Decomposition Temperature</b> | No Data Available |
| <b>Evaporation Rate</b>                               | No Data Available       | <b>Viscosity</b>                 | No Data Available |
| <b>Flammability (Solid, Gas)</b>                      | No Data Available       | <b>Explosive Properties</b>      | No Data Available |
| <b>Upper / Lower Flammability Or Explosive Limits</b> | No Data Available       | <b>Oxidizing Properties</b>      | No Data Available |

**SECTION 10 - STABILITY AND REACTIVITY**

**Stability:** Stable under recommended storage conditions.

**Hazardous Decomposition Products/ Hazardous Polymerization:** No Data Available.

**Incompatibilities:** Strong Oxidizing Agent.

**SECTION 11 - TOXICOLOGICAL INFORMATION**

**Acute Toxicity:**

| <b>Component</b>             | <b>LD50 Oral</b>  | <b>LD50 Dermal</b>  | <b>LD50 Inhalation</b> |
|------------------------------|-------------------|---------------------|------------------------|
| Sodium Dodecyl Sulfate (SDS) | 977 mg/kg (Rat)   | > 2,000 mg/kg (Rat) | 1.51 mg/L (Rat)        |
| Sodium Deoxycholate          | 1,370 mg/kg (Rat) | Not Listed          | Not Listed             |

**SECTION 12 - ECOLOGICAL INFORMATION**



**Ecotoxicity:** Do not empty into drains.

**Persistence and degradability:** Biodegradability: aerobic – Exposure time 28 d. Result: 95% - Readily biodegradable. (OECD Test Guideline 301B)

**Bioaccumulative potential:** Bioaccumulation: Cyprinus carpio (Carp) – 72 h. Bioconcentration factor (BCF): 3.9 – 5.3

**Mobility in soil:** No data available

**Results of PBT and vPvB assessment:** No data available

**Other adverse effects:** No data available

**SECTION 13 - DISPOSAL CONSIDERATIONS**

Whatever cannot be saved for recovery or recycling should be managed in an appropriate and approved waste disposal facility. Processing, use or contamination of this product may change the waste management options. State and local disposal regulations may differ from federal disposal regulations. Dispose of container and unused contents in accordance with federal, state and local requirements.

**SECTION 14 - TRANSPORT INFORMATION**

| DOT           | TDG           | IATA          | IMDG/IMO      |
|---------------|---------------|---------------|---------------|
| Not regulated | Not regulated | Not regulated | Not regulated |

**SECTION 15 - REGULATORY INFORMATION**

**USA Federal Regulation**

**SARA 313/SARA 311/312 HAZARDOUS CATEGORIZATION:** Not Applicable

**Clean Water Act/ Clean Air Act:** No Data Available

**TSCA 12(b)/ OSHA/ CERCLA:** Not Applicable

**California Proposition 65:** This product does not contain any Proposition 65 chemicals.

**STATE RIGHT TO KNOW**

| Massachusetts | New Jersey | Pennsylvania | Illinois | Rhode Island |
|---------------|------------|--------------|----------|--------------|
| -             | Listed     | Listed       | -        | -            |

**HMIS Rating**

| Health hazard | Chronic Health Hazard | Flammability | Physical Hazard |
|---------------|-----------------------|--------------|-----------------|
| 0             |                       | 0            | 0               |

**NFPA Rating**

| Health Hazard | Fire Hazard | Reactivity Hazard |
|---------------|-------------|-------------------|
| 0             | 0           | 0                 |

**US DEPARTMENT OF TRANSPORTATION**

| REPORTABLE QUANTITY (RQ) | DOT MARINE POLLUTANT | DOT SEVER MARINE POLLUTANT |
|--------------------------|----------------------|----------------------------|
| No                       | No                   | No                         |

**U.S. DEPARTMENT OF HOMELAND SECURITY:** This product does not contain any DHS chemicals.

**Canada:** This product has been classified in accordance with the hazard criteria of the Controlled Products Regulations (CPR) and the SDS contains all the information required by the CPR.

**This SDS complies with the requirements of Regulation (EC).**



**Boston BioProducts, Inc**

425 Fortune Blvd.

Milford, MA 01757



**SECTION 16 - OTHER INFORMATION**

The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. Boston BioProducts, Inc. shall not be held liable for any damage resulting from handling or from contact with the above product. This product is sold for laboratory research and development purposes use only.

**Revised: 2023**