



## SAFETY DATA SHEET

### **SECTION 1 – PRODUCT AND COMPANY IDENTIFICATION**

- |            |                                         |                                                                  |
|------------|-----------------------------------------|------------------------------------------------------------------|
| <b>1.1</b> | PRODUCT NAME:<br>CATALOG #:             | PMSF Solution (0.1 M)<br>BP-481                                  |
| <b>1.2</b> | RECOMMENDED USE:                        | Research use or further manufacturing                            |
| <b>1.3</b> | COMPANY IDENTIFICATION:                 | Boston BioProducts, Inc.<br>425 Fortune Blvd., Milford, MA 01757 |
| <b>1.4</b> | EMERGENCY CONTACT:<br>TOLL FREE NUMBER: | 1-800-535-5053<br>1-888-881-8691                                 |

### **SECTION 2 - HAZARDS IDENTIFICATION**

**2.1 Classification of the substance or mixture**  
**GHS Classification in accordance with 29 CFR 1910 (OSHA HCS)**

Flammable liquids	Category 2
Serious eye damage/Eye irritation	Category 2
Specific target organ toxicity (single exposure)	Category 3
Skin Corrosion/Irritation	Category 1B

**2.2 GHS Label elements, including Hazard and Precautionary Statement(s)**  
**Pictogram**



**Signal word: Danger**

**Hazard statement(s)**

H225	Highly flammable liquid and vapor
H314	Causes severe skin burns and eye damage
H319	Causes serious eye irritation.
H336	May cause drowsiness or dizziness.

**Prevention, Response, Storage and Disposal Precautionary Statement(s)**

P210	Keep away from heat/sparks/open flames/hot surfaces. - No smoking
P233	Keep container tightly closed
P240	Ground/bond container and receiving equipment
P260	Do not breathe dust/fume/gas/mist/vapours/spray



P261	Avoid breathing dust/fume/gas/mist/vapors/spray
P264	Wash skin, face, hand thoroughly after handling.
P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.
P301 + P330 + P331	IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
P303 + P361 + P353	IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.
P304 + P340 + P310	IF INHALED: Remove person to fresh air and keep comfortable for breathing. Immediately call a POISON CENTER/doctor.
P305 + P351 + P338 + P310	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor.
P405	Store locked up.
P501	Dispose of contents/ container to an approved waste disposal plant.

### **SECTION 3- COMPOSITION/INFORMATION ON INGREDIENTS**

**Composition:** 100 mM Phenylmethylsulfonyl Fluoride (PMSF) in Isopropyl Alcohol.

<b>Component</b>	<b>CAS Number</b>	<b>Concentration (%)</b>
Phenylmethylsulfonyl Fluoride (PMSF)	329-98-6	1.72%
Isopropyl Alcohol (99% Solution)	67-63-0	Rest

### **SECTION 4 - FIRST-AID MEASURES**

#### **4.1 Description of first aid measures**

**General advice:** Consult a doctor and show this safety data sheet.

- i. If inhaled:** Remove to fresh air and monitor breathing. If breathing becomes difficult, give oxygen. If breathing stops, give artificial respiration. Consult a doctor.
- ii. In case of skin contact:** Immediately wash skin with copious amounts of soap and water for at least 15 minutes. Remove contaminated clothing and shoes and wash before reuse. Consult a doctor.
- iii. In case of eye contact:** Flush with copious amounts of water for at least 15 minutes. Consult a doctor.
- iv. If swallowed:** Rinse mouth with water. Do not induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Consult a doctor.

**4.2 Most important symptoms and effects, both acute and delayed:** To the best of our knowledge, the chemical, physical and toxicological properties have not been thoroughly investigated.

**4.3 Indication of immediate medical attention and special treatment needed:** Show this safety data sheet to the doctor in attendance. Immediate medical attention is required.

### **SECTION 5 - FIRE FIGHTING MEASURES**

#### **5.1 Extinguishing media**

**Suitable extinguishing media:** Use water spray, alcohol-resistant foam, dry chemical, or carbon dioxide

**5.2 Special hazards arising from the substance or mixture:** No Data Available



**5.3 Precautions for fire-fighters:** Wear suitable protective clothing to prevent contact with skin and eyes and self-contained breathing apparatus.

**5.4 Further information:** No data available

## **SECTION 6 - ACCIDENTAL RELEASE MEASURES**

**6.1 Personal precautions, protective equipment, and emergency procedures:** Do not take action without suitable protective clothing - see section 8 of SDS. Evacuate personnel to safe areas. Ensure adequate ventilation. Avoid breathing vapors, mist, dust, or gas.

**6.2 Environmental precautions:** Do not let product enter drains.

**6.3 Methods and materials for containment and cleaning up:** Cover spillage with suitable absorbent material. Using non-spark tools, sweep up material and place in an appropriate container. Decontaminate spill site with 10% caustic solution and ventilate area until after disposal is complete. Hold all material for appropriate disposal as described under section 13 of SDS.

**6.4 Reference to other sections:** For required PPE see section 8. For disposal see section 13.

## **SECTION 7 - HANDLING AND STORAGE**

**7.1 Precautions for safe handling:** Use in a chemical fume hood, with air supplied by an independent system. Avoid inhalation, contact with eyes, skin, and clothing. Avoid the formation of dust and aerosols. Use in a well-ventilated area. Keep away from sources of ignition. Avoid prolonged or repeated exposure.

**7.2 Conditions for safe storage, including any incompatibilities:** Store in cool, well-ventilated area. Keep away from direct sunlight. Keep container tightly sealed until ready for use. Keep in a dry place.

**7.3 Specific end use(s):** Use in a laboratory fume hood where possible

## **SECTION 8 - EXPOSURE CONTROLS/PERSONAL PROTECTION**

### **8.1 Control parameters**

**Components with workplace control parameters:** Contains no substances with occupational exposure limit values.

### **8.2 Exposure Controls**

- i. **Appropriate engineering controls:** Use in a fume hood where applicable. Ensure all engineering measures described under section 7 of SDS are in place. Ensure laboratory is equipped with a safety shower and eye wash station.

### **8.3 Personal protective equipment**

- i. **Eye/face protection:** Face shield and safety glasses Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).
- ii. **Skin protection:** Wear appropriate protective eyeglasses or chemical safety goggles as described by OSHA's eye and face protection regulations in 29 CFR 1910.133 or European Standard EN166
- iii. **Body Protection:** Wear appropriate protective clothing. Complete suit protecting against chemicals, the type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.
- iv. **Respiratory Protection:** Follow the OSHA respirator regulations found in 29 CFR 1910.134 or European Standard EN 149. Use a NIOSH/MSHA or European Standard EN149 approved respirator if exposure limits are exceeded or if irritation or other symptoms are experienced.
- v. **Control of environmental exposure:** Do not let product enter drains.



**SECTION 9 - PHYSICAL AND CHEMICAL PROPERTIES**

**Information on basic physical and chemical properties**

<b>Appearance</b>	Clear, Colorless Liquid	<b>Vapor Pressure</b>	No Data Available
<b>Odor</b>	Alcohol-like	<b>Vapor Density</b>	No Data Available
<b>Odor Threshold</b>	No Data Available	<b>Relative Density</b>	No Data Available
<b>pH</b>	No Data Available	<b>Water Solubility</b>	No Data Available
<b>Melting / Freezing Point</b>	No Data Available	<b>Partition Coefficient</b>	No Data Available
<b>Initial Boiling Point Range</b>	No Data Available	<b>Auto-Ignition Temperature</b>	No Data Available
<b>Flash Point</b>	No Data Available	<b>Decomposition Temperature</b>	No Data Available
<b>Evaporation Rate</b>	No Data Available	<b>Viscosity</b>	No Data Available
<b>Flammability (Solid, Gas)</b>	No Data Available	<b>Explosive Properties</b>	No Data Available
<b>Upper / Lower Flammability or Explosive Limits</b>	No Data Available	<b>Oxidizing Properties</b>	No Data Available

**SECTION 10 - STABILITY AND REACTIVITY**

**Stability:** Stable under recommended storage conditions.

**Hazardous Decomposition Products/ Hazardous Polymerization:** No Data Available.

**Incompatibilities:** Strong Oxidizing Agents.

**SECTION 11 - TOXICOLOGICAL INFORMATION**

**Acute Toxicity:**

<b>Component</b>	<b>LD50 Oral</b>	<b>LD50 Dermal</b>	<b>LC50 Inhalation</b>
Isopropyl alcohol	5045 mg/kg (Rat) 3600 mg/kg (Mouse)	12800 mg/kg (Rat)	72.6 mg/L (Rat) 4 h
Phenylmethylsulfonyl fluoride	200 mg/kg (Rat)	150 mg/kg (Rat)	-

**Skin corrosion/irritation:**

Eyes – Rabbit

Result: Causes severe burns. – 24 h

**Serious eye damage/eye irritation:**

Eyes – Rabbit

Result: Corrosive – 24 h

**Respiratory or skin sensitization:**

Will not occur.

**Germ cell mutagenicity:**

No data available.

**Carcinogenicity:**

IARC	No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible, or confirmed human carcinogen by IARC.
NTP	No component of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.
OSHA	No component of this product present at levels greater than or equal to 0.1% is on OSHA's list of regulated carcinogens.



**Reproductive toxicity:**

No data available.

**Specific target organ toxicity – single exposure:**

No data available.

**Specific target organ toxicity – repeated exposure:**

No data available.

**Aspiration hazard:**

No data available.

**Additional Information:**

RTECS: WB4900000

Material is extremely destructive to tissue of the mucous membranes and upper respiratory tract, eyes, and skin.

**SECTION 12 - ECOLOGICAL INFORMATION**

**Ecotoxicity:**

<b>Component</b>	<b>Freshwater Algae</b>	<b>Freshwater Fish</b>	<b>Microtox</b>	<b>Water Flea</b>
Isopropyl Alcohol	EC50: > 1000 mg/L, 72h (Desmodesmus subspicatus) EC50: > 1000 mg/L, 96h (Desmodesmus subspicatus)	LC50: = 9640 mg/L, 96h flow-through (Pimephales promelas) LC50: > 1400000 µg/L, 96h (Lepomis macrochirus) LC50: = 11130 mg/L, 96h static (Pimephales promelas) LC50: = 10000000 µg/L, 96h (Daphnia)	= 35390 mg/L EC50 Photobacterium phosphoreum 5 min	13299 mg/L EC50 = 48 h9714 mg/L EC50 = 24 h

**Persistence and degradability:**

The methods for determining the biological degradability are not applicable to inorganic substances.

**Bioaccumulative potential/ Mobility in soil:** No data available

**Results of PBT and vPvB assessment:**

PBT/vPvB assessment not available as chemical safety assessment not required/not conducted.

**Other adverse effects:**

An environmental hazard cannot be excluded in the event of unprofessional handling or disposal.

Harmful to aquatic life.



**Chronic Toxicity:** There are no known carcinogenic chemicals in this product.

**SECTION 13 - DISPOSAL CONSIDERATIONS**

Whatever cannot be saved for recovery or recycling should be managed in an appropriate and approved waste disposal facility. Processing, use or contamination of this product may change the waste management options. State and local disposal regulations may differ from federal disposal regulations. Dispose of container and unused contents in accordance with federal, state and local requirements.



**SECTION 14 - TRANSPORT INFORMATION**

DOT	TDG	IATA	IMDG/IMO
<b>UN Number:</b> 1823 <b>Proper Shipping Name:</b> Sodium hydroxide, solid <b>Class:</b> 8 <b>Packing Group:</b> II <b>Reportable Quantity:</b> 1000 lbs <b>Poison Inhalation Hazard:</b> No  	Not regulated	<b>UN Number:</b> 1823 <b>Proper Shipping Name:</b> Sodium hydroxide, solid <b>Class:</b> 8 <b>Packing Group:</b> II  	Not regulated

**SECTION 15 - REGULATORY INFORMATION**

**USA Federal Regulation**

**SARA 302 Components:**

No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

**SARA 313 Components:**

This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

**SARA 311/312 HAZARDOUS CATEGORIZATION:**

No SARA Hazards

**Clean Water Act/ Clean Air Act:** No Data Available

**Clean Air Act:** Not Applicable

**TSCA 12(b)/ OSHA/ CERCLA:** Not Applicable

**California Proposition 65:** This product does not contain any chemicals known to State of California to cause cancer, birth defects, or any other reproductive harm.

**STATE RIGHT TO KNOW**

Component	Massachusetts	New Jersey	Pennsylvania	Illinois	Rhode Island
	Listed	Listed	Listed	-	-

**HMIS Rating**

Health hazard	Chronic Health Hazard	Flammability	Physical Hazard
0		0	0

**NFPA Rating**

Health Hazard	Fire Hazard	Reactivity Hazard
0	0	0



**Boston BioProducts, Inc**

425 Fortune Blvd.

Milford, MA 01757



**US DEPARTMENT OF TRANSPORTATION**

<b>REPORTABLE QUANTITY (RQ)</b>	<b>DOT MARINE POLLUTANT</b>	<b>DOT SEVER MARINE POLLUTANT</b>
<b>No</b>	<b>No</b>	<b>No</b>

**U.S. DEPARTMENT OF HOMELAND SECURITY:** This product does not contain any DHS chemicals.

**Canada:** This product has been classified in accordance with the hazard criteria of the Controlled Products Regulations (CPR) and the SDS contains all the information required by the CPR.

**This SDS complies with the requirements of Regulation (EC).**

**SECTION 16 - OTHER INFORMATION**

The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. Boston BioProducts, Inc shall not be held liable for any damage resulting from handling or from contact with the above product. This product is sold for laboratory research and development purposes use only.

**Revised: 2024**