



## SAFETY DATA SHEET

### **SECTION 1 – PRODUCT AND COMPANY IDENTIFICATION**

- |            |   |   |
|------------|---|---|
| <b>1.1</b> | PRODUCT NAME:<br>CATALOG #:             | Triton X-100<br>P-924   |
| <b>1.2</b> | RECOMMENDED USE:                        | Research Use or Further Manufacturing                           |
| <b>1.3</b> | COMPANY IDENTIFICATION:                 | Boston BioProducts, Inc.<br>425 Fortune Blvd. Milford, MA 01757 |
| <b>1.4</b> | EMERGENCY CONTACT:<br>TOLL FREE NUMBER: | 1-800-535-5053<br>1-888-881-8691                                |

### **SECTION 2 - HAZARDS IDENTIFICATION**

#### **2.1 Classification of the substance or mixture**

#### **GHS Classification in accordance with 29 CFR 1910 (OSHA HCS)**

H302	Harmful if swallowed	Category 4
H315	Causes skin irritation	Category 2
H318	Causes serious eye damage	Category 1
H410	Very toxic to aquatic life with long-lasting effects	Category 1

#### **2.2 GHS Label elements, including Hazard and Precautionary Statement(s)**

##### **Pictogram**



**Signal word:** Danger

##### **Hazard statement(s):**

H302	Harmful if swallowed
H315	May cause skin irritation
H318	Causes serious eye damage
H400	Very toxic to aquatic life
H410	Very toxic to aquatic life with long lasting effects.

##### **Precautionary statements:**

P101	If medical advice is needed, have product container or label at hand
P102	Keep out of reach of children.
P264	Wash skin thoroughly after handling.
P273	Avoid release to the environment.
P280	Wear protective gloves/ eye protection/ face protection.
P301 + P312 + P330	IF SWALLOWED: Call a POISON CENTER/ doctor if you feel unwell. Rinse mouth.



P302 + P352	IF ON SKIN: Wash with plenty of soap and water.
P305 + P351 + P338 + P310	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/ doctor.
P332 + P313	If skin irritation occurs: Get medical advice/ attention.
P362	Take off contaminated clothing and wash before reuse
P501	Dispose of contents/container to an approved waste disposal plant.

**SECTION 3- COMPOSITION/INFORMATION ON INGREDIENTS**

**Composition:** >99% (v/v) Triton X-100 (t-Octylphenoxypolyethoxyethanol 4-(1,1,3,3-Tetramethylbutyl)phenyl polyethylene glycol Polyethylene glycol tert-octylphenyl ether).

Component	CAS Number	Concentration
Triton X-100	9036-19-5	>99% (v/v)

**SECTION 4 - FIRST-AID MEASURES**

**4.1 Description of First Aid Measures**

**General advice:** Consult a doctor and show this safety data sheet.

- i. **If inhaled:** Remove to fresh air and monitor breathing. If breathing becomes difficult, give oxygen. If breathing stops, give artificial respiration. Consult a doctor.
- ii. **In case of skin contact:** Immediately wash skin with copious amounts of soap and water for at least 15 minutes. Remove contaminated clothing and shoes and wash before reuse. Consult a doctor.
- iii. **In case of eye contact:** Flush with copious amounts of water for at least 15 minutes. Consult a doctor.
- iv. **If swallowed:** Rinse mouth with water. Do not induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Consult a doctor.

**4.2 Most important symptoms and effects, both acute and delayed:** To the best of our knowledge, the chemical, physical and toxicological properties have not been thoroughly investigated.

**4.3 Indication of immediate medical attention and special treatment needed:** Show this safety data sheet to the doctor in attendance. Immediate medical attention is required.

**SECTION 5 - FIRE FIGHTING MEASURES**

**5.1 Extinguishing media**

**Suitable extinguishing media:** Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide

**5.2 Special hazards arising from the substance or mixture:** No Data Available.

**5.3 Precautions for fire-fighters:** Wear suitable protective clothing to prevent contact with skin and eyes and self-contained breathing apparatus.

**5.4 Further information:** No data available

**SECTION 6 - ACCIDENTAL RELEASE MEASURES**

**6.1 Personal precautions, protective equipment and emergency procedures:** Do not take action without suitable protective clothing - see section 8 of SDS. Evacuate personnel to safe areas. Ensure adequate ventilation. Avoid breathing vapors, mist, dust or gas.

**6.2 Environmental precautions:** Do not let product enter drains.



**6.3 Methods and materials for containment and cleaning up:** Cover spillage with suitable absorbent material. Using non-spark tools, sweep up material and place in an appropriate container. Decontaminate spill site with 10% caustic solution and ventilate area until after disposal is complete. Hold all material for appropriate disposal as described under section 13 of SDS.

**6.4 Reference to other sections:** For required PPE see section 8. For disposal see section 13.

**SECTION 7 - HANDLING AND STORAGE**

**7.1 Precautions for safe handling:** Use in a chemical fume hood, with air supplied by an independent system. Avoid inhalation, contact with eyes, skin and clothing. Avoid the formation of dust and aerosols. Use in a well-ventilated area. Keep away from sources of ignition. Avoid prolonged or repeated exposure.

**7.2 Conditions for safe storage, including any incompatibilities:** Store in cool, well-ventilated area. Keep away from direct sunlight. Keep container tightly sealed until ready for use. Keep in a dry place.

**7.3 Specific end use(s):** Use in a laboratory fume hood where possible

**SECTION 8 - EXPOSURE CONTROLS/PERSONAL PROTECTION**

**8.1 Control parameters**

**Components with workplace control parameters:** Contains no substances with occupational exposure limit values.

**8.2 Exposure Controls**

- i. **Appropriate engineering controls:** Use in a fume hood where applicable. Ensure all engineering measures described under section 7 of SDS are in place. Ensure laboratory is equipped with a safety shower and eye wash station.

**8.3 Personal protective equipment**

- i. **Eye/face protection:** Face shield and safety glasses Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).
- ii. **Skin protection:** Wear appropriate protective eyeglasses or chemical safety goggles as described by OSHA's eye and face protection regulations in 29 CFR 1910.133 or European Standard EN166
- iii. **Body Protection:** Wear appropriate protective clothing. Complete suit protecting against chemicals, the type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.
- iv. **Respiratory Protection:** Follow the OSHA respirator regulations found in 29 CFR 1910.134 or European Standard EN 149. Use a NIOSH/MSHA or European Standard EN149 approved respirator if exposure limits are exceeded or if irritation or other symptoms are experienced.
- v. **Control of environmental exposure:** Do not let product enter drains

**SECTION 9 - PHYSICAL AND CHEMICAL PROPERTIES**

**Information on basic physical and chemical properties**

<b>Appearance</b>	Clear, Colorless Liquid	<b>Vapor Pressure</b>	No Data Available
<b>Odor</b>	No Data Available	<b>Vapor Density</b>	No Data Available
<b>Odor Threshold</b>	No Data Available	<b>Relative Density</b>	No Data Available
<b>pH</b>	No Data Available	<b>Water Solubility</b>	No Data Available
<b>Melting/Freezing Point</b>	No Data Available	<b>Partition Coefficient</b>	No Data Available
<b>Initial Boiling Point Range</b>	No Data Available	<b>Auto-Ignition Temperature</b>	No Data Available
<b>Flash Point</b>	No Data Available	<b>Decomposition Temperature</b>	No Data Available
<b>Evaporation Rate</b>	No Data Available	<b>Viscosity</b>	No Data Available
<b>Flammability (Solid, Gas)</b>	No Data Available	<b>Explosive Properties</b>	No Data Available



<b>Upper/Lower Flammability or Explosive Limits</b>	No Data Available	<b>Oxidizing Properties</b>	No Data Available
---	-------------------	-----------------------------	-------------------

**SECTION 10 - STABILITY AND REACTIVITY**

**Stability:** Stable under recommended storage conditions.

**Hazardous Decomposition Products/ Hazardous Polymerization:** No data available

**Incompatibilities:** Strong oxidizing agents.

**SECTION 11 - TOXICOLOGICAL INFORMATION**

**Acute Toxicity:**

<b>Acute toxicity</b>	<b>Skin corrosion/irritation</b>	<b>Serious eye damage/eye irritation</b>	<b>Respiratory or skin sensitization</b>
LD50 Oral - Rat - 1,900-5,000 mg/kg LD50 Dermal - Rabbit - >3,000 mg/kg	Skin - Rabbit Result: irritating - 4 h (OECD Test Guideline 404) Remarks: The value is given in analogy to the following substances: 4-(1,1,3,3-tetramethylbutyl)phenol	Eyes - Rabbit Result: Risk of serious damage to eyes. (Draize Test) Remarks: Risk of corneal clouding.	Sensitisation test: - Human Result: negative Remarks: (External MSDS) Patch test on human volunteers did not demonstrate sensitization properties

**SECTION 12 - ECOLOGICAL INFORMATION**

**Ecotoxicity:**

<b>Toxicity to fish</b>	<b>Toxicity to daphnia and other aquatic invertebrates</b>	<b>Toxicity to algae</b>	<b>Toxicity to fish (chronic toxicity)</b>	<b>Toxicity to daphnia and other aquatic invertebrates (Chronic toxicity)</b>
LC50 - Pimephales promelas (fathead minnow) - 4 - 8.9 mg/l - 96 h  semi-static test LC50 - Leuciscus idus (Golden orfe) - 0.26 mg/l - 96 h (OECD Test Guideline 203)  Remarks: The value is given in analogy to the following substances: 4-(1,1,3,3-tetramethylbutyl)phenol	LC50 - Daphnia magna (Water flea) - 18 - 26 mg/l - 48 h  static test EC50 - Daphnia magna (Water flea) - 0.011 mg/l - 48 h Remarks: (ECOTOX Database) The value is given in analogy to the following substances: 4-(1,1,3,3-tetramethylbutyl)phenol	Growth inhibition LOEC - Pseudokirchneriella subcapitata - 16.0 mg/l - 96 h Growth inhibition NOEC - Pseudokirchneriella subcapitata - 8.0 mg/l - 96 h	flow-through test - Danio rerio (zebra fish) - 0.012 mg/l (OECD Test Guideline 210) Remarks: The value is given in analogy to the following substances: 4-(1,1,3,3-tetramethylbutyl)phenol	semi-static test NOEC - Daphnia magna (Water flea) - 0.03 mg/l - 21 d (OECD Test Guideline 202) Remarks: The value is given in analogy to the following substances: 4-(1,1,3,3-tetramethylbutyl)phenol

**Persistence and degradability:** aerobic - Exposure time 28 d Result: 22 % - Not readily biodegradable. (OECD Test Guideline 301C)

**Bioaccumulative potential/ Mobility in soil:**

**Results of PBT and vPvB assessment:** No data available

**Other adverse effects:** No data available

**Chronic Toxicity:** There are no known carcinogenic chemicals in this product



**SECTION 13 - DISPOSAL CONSIDERATIONS**

Whatever cannot be saved for recovery or recycling should be managed in an appropriate and approved waste disposal facility. Processing, use or contamination of this product may change the waste management options. State and local disposal regulations may differ from federal disposal regulations. Dispose of container and unused contents in accordance with federal, state and local requirements.

**SECTION 14 - TRANSPORT INFORMATION**

<b>DOT</b>	<b>TDG</b>	<b>IATA</b>	<b>IMDG/IMO</b>
<b>UN#: 3082</b> <b>Proper Shipping Name:</b> ENVIRONMENTALLY HAZARDOUS SUBSTANCE, LIQUID, N.O.S. (Triton X-100) <b>Hazard Class: 9</b> <b>Packaging Group: III</b>	<b>UN#: 3082</b> <b>Proper Shipping Name:</b> ENVIRONMENTALLY HAZARDOUS SUBSTANCE, LIQUID, N.O.S. (Triton X-100) <b>Hazard Class: 9</b> <b>Packaging Group: III</b>	<b>UN#: 3082</b> <b>Proper Shipping Name:</b> ENVIRONMENTALLY HAZARDOUS SUBSTANCE, LIQUID, N.O.S. (Triton X-100) <b>Hazard Class: 9</b> <b>Packaging Group: III</b>	<b>UN#: 3082</b> <b>Proper Shipping Name:</b> ENVIRONMENTALLY HAZARDOUS SUBSTANCE, LIQUID, N.O.S. (Triton X-100) <b>Hazard Class: 9</b> <b>Packaging Group: III</b>

**SECTION 15 - REGULATORY INFORMATION**

**USA Federal Regulation**

**SARA 313:** Not Applicable

**SARA 311/312 HAZARDOUS CATEGORIZATION:**

<b>Acute Health Hazard</b>	<b>Chronic Health Hazard</b>	<b>Fire Hazard</b>	<b>Sudden Release of Pressure Hazard</b>	<b>Reactive Hazard</b>
Yes	No	No	No	No

**Clean Water Act/ Clean Air Act:** No Data Available

**Clean Air Act:** Not Applicable

**TSCA 12(b)/ OSHA/ CERCLA:** Not Applicable

**California Proposition 65:** This product does not contain any Proposition 65 chemicals.

**STATE RIGHT TO KNOW**

<b>Massachusetts</b>	<b>New Jersey</b>	<b>Pennsylvania</b>	<b>Illinois</b>	<b>Rhode Island</b>
-	Listed	Listed	-	-

**HMIS Rating**

<b>Health hazard</b>	<b>Chronic Health Hazard</b>	<b>Flammability</b>	<b>Physical Hazard</b>
2		1	2

**NFPA Rating**

<b>Health Hazard</b>	<b>Fire Hazard</b>	<b>Reactivity Hazard</b>
2	1	0

**US DEPARTMENT OF TRANSPORTATION**

<b>REPORTABLE QUANTITY (RQ)</b>	<b>DOT MARINE POLLUTANT</b>	<b>DOT SEVER MARINE POLLUTANT</b>
No	No	No



**Boston BioProducts, Inc**

425 Fortune Blvd.

Milford, MA 01757

**U.S. DEPARTMENT OF HOMELAND SECURITY:** This product does not contain any DHS chemicals.

**Canada:** This product has been classified in accordance with the hazard criteria of the Controlled Products Regulations (CPR) and the SDS contains all the information required by the CPR.

**This SDS complies with the requirements of Regulation (EC).**

#### **SECTION 16 - OTHER INFORMATION**

The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. Boston BioProducts, Inc shall not be held liable for any damage resulting from handling or from contact with the above product. This product is sold for laboratory research and development purposes use only.

**Revised: 2024**